

May 30th, 2023

2023\_ILT\_04

Specimen	Status	Location	Conc	Date	Notes	Volume	Lym	Lym+Mon	Total
2000cc- 65 Inf070-7 a-5	HU ♀ NWCC	Box 3 AC1	~5 45.37 (5)			0.75	12.8	18.0	9.6
1350cc- 65 a-6	RWCC		46.87 (9)		open, 102 inside, Haw issue	0.75	10.6	15.5	7.95
1800cc- 65 Inf374-4 a-9	HEU-lo SV50	Box 9B H8	~7.5 46.67 (1)			0.75	9.5	13.6	7.12
1780cc- 65 a-10	SV51		46.17 (7)			0.75	8.1	12.6	6.08
65 Inf287-7 a-4	HEU-lo ♀ SLZ1	Box 7B B10	~5 (3)			0.75	4.8		4.70
1460cc- 65 a-7	SLZ2		46.17 (8)		100% tape! Horvitz	0.75	6.3	11.7	4.70
1450cc- 65 Inf191-1 a-5	HEU-hi ♀ SLAN	Box 6B C7	~10 (6)		Vapour, worse led?	1	7.7	13.2	
800cc- 40 ND050	Adult		15EC 41.47 (2)	1/13/23		1	8.5	12.2	
830cc- 40 ND006	Norm		15EC 48.27 (1)	1/20/23		1	11.8	14.7	

10:18 am Adult Trans 10:41 DMSO adult  
 10:31 am 10:47 spin 11:04 DMSO  
 10:51 am 11:04 DMSO  
 11:39 am stain for count  
 12:09 pm count done  
 12:29 aliquots  
 1:08 pm R10's  
 1:24 ND006

n=17      n=10

0000	00
0000	00
000	00
2.5M 000	00
000	00

1:43pm Cells into Incubator

→ 7:43pm

Note Inf191 vs PMA ↑ pma dose  
 - vs ND050 vs PMA  
 maybe up to x2? Not missing  
 @ this point, but note if after rank  
 ↑ 1.44 raised 830/1 added 15/1000  
 to restore RNA/ml balance, at cost  
 cells/PL conc → 1.5M

5:38 most Abs & reagents prepped  
 Exhaustive search for BV420 CO161, CO2C, CO141  
 6:03 everything - zombie W2 prepped //

7:51 pm spin @

H/O @ 8:00 - spin @ 8:29 pm

20:24 sc spin

Hot samples @ 8:47 → 9:17 pm Spin @ 9:21

n=36 sc's

Zombie sc's @ 8:52 → 9:07 pm  
 hot sc's @ 8:59 → 9:29  
 Cold sc's @ 9:08 → 9:38  
 9:32 tets prepped

Tets @ 9:48 → 10:28 pm  
 Abs @ 10:00 → 10:30 pm

Sc Fix Perm @ 9:56 → 10:06 → 10:16

sc first spin @ 10:22

Cold samples @ 22:49 → 23:29

sc's intrac @ 10:58 → 11:38 pm

sc's done @ 11:06 pm

Rbc lysis @ 11:26 pm

↓

Spin @ 11:34 pm

Mark specimens

Sample FoxPerm @ 11:48 → 58 → 00:08

SES done @ 11:50 pm w/ 20  $\mu$ l of 0.4% PFA

→ 1st wash @ 12:12 am

2nd wash @ 12:24 am

Intra @ 12:37 am → 1:17 am


1:19 am final spin

Done @ 1:33 am, placed in 30  $\mu$ l

In cell ventral

is  $\sigma$  of a mother expression?


Cobran in the lines



NO50 35  $\mu$ l  
37  $\mu$ l

Start @ 8:30 pm < 35  $\mu$ l > actually a ctrl well this time.  
200,000


9:45 instants



Inf 070 ctrl

4900/92 55  $\mu$ l

PMD 4100/68 45  $\mu$ l



Inf 374

6200/147 45  $\mu$ l  
6100/140 43  $\mu$ l



Inf 287 43  $\mu$ l  
(poorer shape) 50  $\mu$ l  
↑ evidence extra AF



Inf 191  
35  $\mu$ l 50  $\mu$ l  
41  $\mu$ l  
Very few alive  
in PMD

10:07 15000/743 48  $\mu$ l  
60  $\mu$ l

2.5 m cells / tube  
events

\* How communicate effectively to peers, SRLs, etc not in same frame/thought? \*

\* Inf 374 82  $\mu$ l tet → hiccups clump 7  $\mu$ l. \*

287 25  $\mu$ l / 3900

\* How much variation paracalling is acceptable? \*

Done @ 11:00 pm

May 30th, 2023

2023\_ILT\_04







49 21 10  
40 24 10

KCR out O

#	Filter	Single color (u)	Ref cell	Unbinding ctrl	Fluorochrome	Marker	Clone	Val lot #	During stimuli	17	U/D 15 min (RT)	Tetramer 40 min @ RT	HisStain 30min @ 37C	17	ColdStain 30min @ 4C	17	RBC Lys. FixPerm	Spiked 40 min @ RT	17
1	UV7				BUV385	CD45L	SK11								1.2	20.4			
2	UV7				AF	AF-LIVE													
3	UV7				BUV486	CD8	RP418								0.7	11.9			
4	UV10				BUV661	CD69	FM50								0.5	8.5		1	17.0
5	UV11				BUV615	CCR4	1C1												
6	UV4				BUV661	CCR4	B6								0.7	11.9			
7	UV16				BUV37	CCR3	1C6/CCR3								1.3	22.1			
8	UV16				BUV605	CD4	SK3												
9	UV16				BUV421	CD127	A01905												
10	UV16				Pacific Blue	CD14	MS2								2.0	34.0		0.5	8.5
11	UV16				BUV480	CD161	SUZC1								2.0	34.0			
12	UV16				BUV310	CD38A	HL100												
13	UV16				BUV605	CD56	5.H111								0.7	11.9			
14	UV16				BUV605	CD57	6043H7											0.1	1.7
15	UV16				BUV711	CD7	M1701											0.7	11.9
16	UV16				BUV750	IRF7	B27												
17	B3				BUV750	CD86	11A9											1.5	25.5
18	B3				ETC/AF488	CD135	8B11												
19	B3				Spark Blue 550	CD3	SK7												
20	B4				PE-Cy5.5	CD26	B55b												
21	B6				PE	NG2D	1011												
22	B8				PE-DiIc488	TNF4	MA611												
23	B11				PE-Cy5	CD25	MA251												
24	B11				PE-CD70	PD1	PD1.1.3												
25	B2				APC	CD16	3C8												
26	B4				Alexa488/AF488	CD107a	3C10												
27	B6				APC-CD70	CD107a	H4A3												
28	B7				Zombie NIR	L/D	N/A												
29	B8				APC/Fire 750	CD27	O323												
					APC/Fire 810	CD38	HL12												
And UNSTAINED CONTROLS III																			
Antibody Total									6.0	102									
R10 Media									14.5	246.5									
Pipette draw volume / sample									19.5										
Antibody Total																			
Brilliant Stain																			
Pipette draw volume / sample																			
Antibody Total									17.8	303	12.5	212.5						9.5	161.5
Brilliant Stain									50	850	50.0	850						11	187
Pipette draw volume / sample									65		59.5							19.5	

UV Mix	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Tetramer Mix	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Tetramer Control Mix	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Zombie NIR	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Tetramer Mix	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Tetramer Control Mix	W22/1418	1.5	10
AF488	W22/1418	1.2	15
AF647	W22/1418	1.2	12
Pipette draw volume / sample		1.78	178

Number Samples	17	136	5.1	5.1	68	0.755
Number Unstained	10	30	3	6	20	0.25
Number Surface StCs (34)	30	60	0	12	24	0.25
Number Intracellular StCs	6	12	0	12	24	0.25
Total	288	8.1	14.3	172	1.915	

LRG	144	144	144
PBS	1.7235	1.7235	1.7235
PFA	0.1515	0.1515	0.1515
RBC	0.81	0.81	0.81
Water	7.29	7.29	7.29

## Simplified Protocol

Thaw cells, DMSO, count.  
Collect, count, aliquot cells 3.0-6.0 Cells R10 / 5ml polypropylene tube  
Bring volume up to 100 µl R10, and 100 µl PFA/CD107a  
Cap and incubate at 37°C for 6 hours

Wash with 2 ml PBS, spin down 1300rpm 8min

800 µl of Lysedbed mix (1:2500) @ RT for 15min

Wash 2 ml 5% PBS-FBS, spin 1300 rpm, 8min

Add HisStain mix, incubate @ 37°C for 30 min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min

Add Tetramers, incubate @ RT for 10 min

Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min

Add ColdStain mix, incubate @ 4°C for 30min

Add 300-500 µl 1x RBC Lysis for 3 minutes

Wash 2 ml 5% PBS-FBS 1400 rpm, 6 min

300 µl BD FixPerm, incubate @ 4°C for 20min

(vortex every 10 minutes)

First PermWash: 1 ml PermWash 1500 rpm 6 min

Second Perm Wash: 1 ml PermWash 1500 rpm 6 min

Add Intracellular stain, incubate @ RT for 40min

First PermWash: 2 ml PermWash 1500 rpm 6 min

Resuspend in 70 µl 0.4% PFA-PBS

Cap tubes, wrap rack in foil, store at 4°C





Spectrum	UV		Violet		Blue		Yellow-Green		Red	
373	UV1									
388	UV2	BUV395	CD62L							
428	UV3			BV421	CD127					
443	UV4									
458	UV5			PacBlue	CD14/19					
473	UV6	AF								
508										
514	UV7	BUV496	CD8	BV480	CD161	B1				
525										
542	UV8					B2				
582	UV9	BUV563	CD69	BV510	FITC/AF488	B3				
598				BV570	SparkBlue 550	B4				
613	UV10	BUV615	CCR4			B5	PE	NKG2D		
664	UV11	BUV661	V62	BV605	CD56	B6	PE-Dazzle594	TNF $\alpha$	APC	CD16
679				BV650	CCR7	B7		CD25	AF647	Va7.2/InMR1
697	UV12					B8				
717	UV13			BV711	CD7	B9	PE-Cy5			
738						B10				
750	UV14	BUV737	CXCR3	BV750	IFN $\gamma$	B11			APC-R700	CD107a
760						B12				
783	UV15					B13			Zombie NIR	Viability
812	UV16	BUV805	CD4	BV786	CCR6	B14	Pe-Vio770	PD1	APC-Fire 750	CD27
									APC-Fire 810	CD38

Spectrum	UV		Violet		Blue		Yellow-Green		Red	
373	UV1									
388	UV2	BUV395	CD62L							
428	UV3			BV421	CD127					
443	UV4									
458	UV5			PacBlue	CD14/19					
473	UV6	AF								
508										
514	UV7	BUV496	CD8	BV480	CD161	B1				
525										
542	UV8					B2				
582	UV9	BUV563	CD69	BV510	FITC/AF488	B3				
598				BV570	SparkBlue 550	B4				
613	UV10	BUV615	CCR4			B5				
664	UV11	BUV661	V62	BV605	CD56	B6	PE-Dazzle594	TNF $\alpha$	APC	CD16
679				BV650	CCR7	B7		CD25	AF647	Va7.2/InMR1
697	UV12					B8				
717	UV13			BV711	CD7	B9	PE-Cy5			
738						B10				
750	UV14	BUV737	CXCR3	BV750	IFN $\gamma$	B11			APC-R700	CD107a
760						B12				
783	UV15					B13			Zombie NIR	Viability
812	UV16	BUV805	CD4	BV786	CCR6	B14	Pe-Vio770	PD1	APC-Fire 750	CD27
									APC-Fire 810	CD38



May 30th, 2023

2023\_ILT\_04





May 30th, 2023

2023\_ILT\_04



May 30th, 2023

2023\_ILT\_04

